


Dateiname / Filename 170801_ChHu_TecDoc_Fibox 4 FTC-PStx_dv1.docx	Specification Sheet	
Created 01.08.2017	Fibox 4 FTC-PStx	Seite / Page 1/2

Technical Documentation

Fibox 4 FTC-PStx

1. Overview


The Fibox 4 FTC-PStx is a completely stand-alone, portable fiber optic oxygen meter with integrated flow-through cell. The 6 mm steel tube of the integrated flow-through measurement cell can easily be connected to perfused systems and different sampling ports of the process line. Gases and liquids can be pumped through the Fibox 4 FTC-PStx enabling easy control of oxygen at different stages of production processes. Oxygen measurements are temperature compensated using an integrated Pt100 temperature sensor. The oxygen meter with integrated flow-through cell is available for wide, trace and ultra-trace measurement ranges. Fibox 4 FTC-PStx is delivered pre-calibrated so measurements can be started right away. This ensures a worry-free workflow while conducting most accurate oxygen measurements.



Device	Measurement Range
Fibox 4 FTC-PSt3	Wide range: 0 – 100 % O ₂ , detection limit 15 ppb
Fibox 4 FTC-PSt6	Trace range: 0 – 5 % O ₂ , detection limit 1 ppb
Fibox 4 FTC-PSt9	Ultra-trace range: 0 – 200 ppmv gaseous O ₂ , detection limit 0.5 ppm

Features:

- Pre-calibrated, ready-to-use
- Temperature compensation with integrated Pt100 temperature sensor
- Measures dissolved and gaseous O₂ (PSt9 gaseous only)
- Portable, stand-alone device with display, battery and data logger
- Splash proof, robust housing
- Compatible with PreSens Measurement Studio 2 (PMS2)

Dateiname / Filename 170801_ChHu_TecDoc_Fibox 4 FTC-PStx_dv1.docx	Specification Sheet	
Created 01.08.2017	Fibox 4 FTC-PStx	Seite / Page 2/2

2. Sensor Characteristics

	PSt3		PSt6		PSt9
Specifications	Gaseous & Dissolved O ₂	Dissolved O ₂	Gaseous & Dissolved O ₂	Dissolved O ₂	Gaseous O ₂
Measurement range	0 – 100 % O ₂ 0 – 1000 hPa	0 – 45 mg/L 0 – 1400 µmol/L	0 – 4.2 % O ₂ 0 – 41.4 µmol/L	0 – 1.8 mg/L 0 – 56.9 µmol/L	0 – 200 ppm
Limit of detection	0.03 % O ₂	15 ppb	0.002 % O ₂	1 ppb	0.5 ppm
Resolution	± 0.01 % O ₂ at 0.21 % O ₂ ± 0.1 % O ₂ at 20.9 % O ₂ ± 0.1 hPa at 2 hPa ± 1 hPa at 207 hPa	± 0.14 µmol/L at 2.83 µmol/L ± 1.4 µmol/L at 238.1 µmol/L	± 0.0007 % O ₂ at 0.002 % O ₂ ± 0.0015 % O ₂ at 0.2 % O ₂ ± 0.007 hPa at 0.023 hPa ± 0.015 hPa at 2.0 hPa	± 0.010 µmol/L at 0.03 µmol/L ± 0.020 µmol/L at 2.8 µmol/L	10 ± 0.5 ppm 100 ± 0.8 ppm 200 ± 1.5 ppm
Accuracy	± 0.4 % O ₂ at 20.9 % O ₂ ± 0.05 % O ₂ at 0.2 % O ₂		± 1 ppb or ± 3 %; whichever is higher		± 2 ppm or ± 5 %; whichever is higher
Measurement temperature range	0 to +50 °C		0 to +50 °C		0 to +40 °C
Response time (t ₉₀)	< 6 sec.	< 40 sec.	< 6 sec.	< 40 sec.	< 3 sec.
Properties					
Compatibility	Aqueous solutions, ethanol, methanol				Gas phase only
No cross-sensitivity with	pH 1 – 14 CO ₂ , H ₂ S, SO ₂ Ionic species				CO ₂ , SO ₂
Cross-sensitivity	Organic solvents, such as acetone, toluene, chloroform or methylene chloride Chlorine gas				Organic vapor, Chlorine gas
Sterilization procedure	Ethylene oxide (EtO)				-
Cleaning procedure	3 % H ₂ O ₂ Acidic agents (HCl, H ₂ SO ₄), max. 4 – 5 % at room temperature				-
Calibration	Two-point calibration with oxygen-free environment (nitrogen, sodium sulfite) and air-saturated environment		Two-point calibration in oxygen-free environment (nitrogen) and a second calibration value optimally between 1 and 2 % oxygen		Two-point calibration in oxygen-free environment (nitrogen 6.0) and a second calibration value optimally between 100 and 200 ppm gaseous oxygen
Storage Stability	2 years provided the sensor material is stored in the dark (-10 to +60 °C)				